



Recombinant Cytolethal distending toxin subunit B (cdtB)

Product Code	CSB-MP671598ENL
Abbreviation	cdtB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46669
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	DL TDFRVATWNL QGASATTESK WNINVRQLIS GENAVDILAV QEAGSPPSTA VDTGTLIPSP GIPVRELIWN LSTNSRPQQV YIYFSAVDAL GGRVNLALVS NRRADEVFVL SPVRQGGRRPL LGIRIGNDAF FTAHAIAMRN NDAPALVEEV YNFFRDSRDP VHQALNWMIL GDFNREPADL EMNLTVPVRR ASEIISPAAA TQTSQRDLDY AVAGNSVAFR PSPLQAGIVY GARRTQISSD HFPVGVSR
Source	Mammalian cell
Target Names	cdtB
Protein Names	Recommended name: Cytolethal distending toxin subunit B Short name= CDT B Alternative name(s): Deoxyribonuclease CdtB EC= 3.1.-.-
Expression Region	19-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.