



Recombinant Escherichia coli Protein SsnA (ssnA)

Product Code	CSB-EP671605ENV
Abbreviation	ssnA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46812
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MLILKNVTAV QLHPAKVQEG VDIAIENDVI VAIGDALTQR YPDASFKEMH GRIVMPGIVC SHNHFYSGLS RGIMANIAPC PDFISTLKNL WWRLDRALDE ESLYYSGLIC SLEAIKSGCT SVIDHHASPA YIGGSLSTLR DAFLKVGLRA MTCFETTDNRN NGIKELQEGV EENIRFARLI DEAKKATSEP YLVEAHIGAH APFTVPDAGL EMLREAVKAT GRGLHIHAAE DLYDVSYSHH WYGKDLLARL AQFDLIDSKT LVAHGLYLSK DDITLLNQRD AFLVHNARSN MNNHVGYNHH LSDIRNLALG TDGIGSDMFE EMKFAFFKHR DAGGPLWPDS FAKALTNGNE LMSRNFGAKF GLLEAGYKAD LTICDYNSPT PLLADNIAGH IAFGMGSGSV HSVMVNGVMV YEDRQFNFDK DSIYAQARKA AASMWRMDA LA
Source	E.coli
Target Names	ssnA
Protein Names	Recommended name: Protein SsnA
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.