



Recombinant Type I restriction enzyme Ecoprrl M protein (hsdM)

Product Code	CSB-EP671725ENL-B
Abbreviation	hsdM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q47163
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MKMTSIQQRA ELHRQIWQIA NDVRGSVDGW DFKQYVLGAL FYRFISENFS SYIEAGDDSI CYAKLDDSVI TDDIKDDAIK TKGYFIYPSQ LFCNVAAKAN TNDRLNADLN SIFVAIESSA YGYPSEADIK GLFADFDTTS NRLGNTVKDK NARLAAVLKG VEGLKLGDFN EHQIDLFGDA YEFLISNYAA NAGKSGGEFF TPQHVSKLIA QLAMHGQTHV NKIYDPAAGS GSLLLQAKKQ FDDHIIIEGF FGQEINHHTY NLARMNMLFH NINYDKFDIK LGNTLTEPHF RDEKPFDAIV SNPPYSVKWI GSDDPTLIND ERFAPAGVLA PKSKADFAFV LHALNYLSAK GRAAIVCFPG IFYRGGAEQK IRQYLVDNNY VETVISLAPN LFFGTTIAVN ILVLSKHKTD TKVQFIDASE LFKKETNNNI LTDAHIEQIM QVFASKEDVA HLAKSVAFET VVANDYNLSV SSYVEAKDTR EIIDIAELNA ELKTTVSKID QLRKDIDAIV AEIEGCEVQK
Source	E.coli
Target Names	hsdM
Protein Names	Recommended name: Type I restriction enzyme Ecoprrl M protein Short name= M.Ecoprrl EC= 2.1.1.72
Expression Region	1-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.