



Recombinant Hylobates lar Glutathione peroxidase 2 (GPX2)

Product Code	CSB-EP672241HXW-B
Abbreviation	GPX2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4AEH9
Product Type	Recombinant Protein
Immunogen Species	Hylobates lar (Common gibbon) (White-handed gibbon)
Purity	>85% (SDS-PAGE)
Sequence	MAFIKSFYD LSAISLDGEK VDFNTFRGRA VLIENVASLU GTTTRDFTQL NELQCRFPRR LVVLGFPCNQ FGHQENCQNE EILNSLKYVR PGGGYQPTFT LVQKCEVNGQ NEHPVFAYLK DKLPYPYDDP FSLMTPKLI IWSPVRRSDV AWNFEKFLIG PEGEPFRRYS RTFPTINIEP DIKRLKVAI
Source	E.coli
Target Names	GPX2
Protein Names	Recommended name: Glutathione peroxidase 2 Short name= GPx-2 Short name= GSHPx-2 EC= 1.11.1.9 Alternative name(s): Glutathione peroxidase-gastrointestinal Short name= GPx-GI Short name= GSHPx-GI
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.