



# Recombinant Didelphis albiventris Fibroblast growth factor 9 (FGF9)

<b>Product Code</b>	CSB-MP672496DBAG
<b>Abbreviation</b>	FGF9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3ZFI5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Didelphis albiventris (White-eared opossum)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LGEVGN Y FGVQDAAVPF GNVPVLPVDS PVLLSDHLGQ SEAGGLPRGP AVTDL D H L K G I L R R R Q L Y C R T G F H L E I F P N G T I Q G T R Q D H S R F G I L E F I S I A V G L V S I R G V D S G L Y L G M N E K G E L Y G S E K L T Q E C V F R E Q F E E N W Y N T Y S S N L Y K H V D T G R R Y Y V A L N K D G T P R E G T R T K R H Q K F T H F L P R P V D L E K V P E L Y K D I L S Q S
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FGF9
<b>Protein Names</b>	Recommended name: Fibroblast growth factor 9 Short name= FGF-9 Alternative name(s): Glia-activating factor Short name= GAF Heparin-binding growth factor 9 Short name= HBGF-9
<b>Expression Region</b>	4-209
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.