



Recombinant Arabidopsis thaliana 60S acidic ribosomal protein P0-2 (RPP0B)

Product Code	CSB-EP672581DOA
Abbreviation	RPP0B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q42112
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVKATKAEKK IAYDTKLCQL IDEYTDQILVV AADNVGSTQL QNIRKGLRGD SVVLMGKNTM MKRSVRIHSE NTGNTAILNL LPLLQGNVGL IFTKGDLKEV SEEVAKYKVG APARVGLVAP IDVVVQPGNT GLDPSQTSFF QVLNIPTKIN KGTVEIITPV ELIKQGDKVG SSEAALLAKL GIRPFSYGLV VQSVYDNGSV FSPEVLDLTDL DQLVEKFASG ISMVTSLALA VSYPTLAAAP HMFINAYKNA LAIIVATEYT FPQAEKVKEY LKDPSKFAVA SVAAVSADAG GGAPAAAKVE EKEESDEEDY GGDFGLFDEE
Source	E.coli
Target Names	RPP0B
Protein Names	Recommended name: 60S acidic ribosomal protein P0-2
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.