



# Recombinant Arabidopsis thaliana Glutamate-1-semialdehyde 2,1-aminomutase 2, chloroplastic (GSA2)

<b>Product Code</b>	CSB-MP672588DOA
<b>Abbreviation</b>	GSA2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q42522
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VSVE EKTKKFTLQK SEEFNAAKN LMPGGVNSPV RAFKSVGGQP VVMDSAKGSR IRDIDGNEYI DYVGSWGPAL IGHADDEVLA ALAETMKKGT SFGAPCLEN VLAEMVISAV PSIEMVRVFN SGTEACMGVL RLARAFTGKQ KFIKFEHCYH GHANSFLVKA GSGVATLGLP DSPGVPKAAT SDTLTAPYND IAAVEKLFEA NKGEIAAIL EPVVGNSGFI TPKPEFIEGI RRITKDNAL LIFDEVMTGF RLAYGGAQEY FGITPDLTTL GKIIGGGLPV GAYGGRRDIM EMVAPAGPMY QAGTLSGNPL AMTAGIHTLK RLSQPGTYEY LDKITKELTN GILEAGKKTG HAMCGGYISG MFGFFFTEGP VYDFSDAKKS DTEKFGKFFR GMLEEGVYLA PSQFEAGFTS LAHTSEDIQF TIAAAEKVLS RL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	GSA2
<b>Protein Names</b>	Recommended name: Glutamate-1-semialdehyde 2,1-aminomutase 2, chloroplastic Short name= GSA 2 EC= 5.4.3.8 Alternative name(s): Glutamate-1-semialdehyde aminotransferase 2 Short name= GSA-AT 2
<b>Expression Region</b>	37-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.