



Recombinant Bordetella pertussis Type III pantothenate kinase (coaX)

Product Code	CSB-MP672674BUA
Abbreviation	coaX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q45338
Product Type	Recombinant Protein
Immunogen Species	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
Purity	>85% (SDS-PAGE)
Sequence	MIILIDSGNS RLKVGWFDPD APQAAREPAP VAFDNLDLDA LGRWLATLPR RPQRALGVNV AGLARGEAIA ATRLRAGGCDI RWLRAQPLAM GLRNGYRNP QLGADRWACM VGVLARQPSV HPPLLVASFG TATTLDTIGP DNVFPGGLIL PGPAMMRGAL AYGTAHLPLA DGLVADYPID THQAIASGIA AAQAGAIVRQ WLAGRQRYGQ APEIYVAGGG WPEVRQEAEER LLAVTGAAFG ATPQPTYLDS PVL DGLAALA A QGAPTA
Source	Mammalian cell
Target Names	coaX
Protein Names	Recommended name: Type III pantothenate kinase EC= 2.7.1.33 Alternative name(s): Bvg accessory factor PanK-III Pantothenic acid kinase
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.