



Recombinant Escherichia coli Uncharacterized protein ygbK (ygbK)

Product Code	CSB-BP672733ENV
Abbreviation	ygbK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46889
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIKIGVIADD FTGATDIASF LVENGLPTVQ INGVPTGKMP EAIDALVISL KTRSCPVEA TQQSLAALSW LQQQGCKQIY FKYCSTFDST AKGNIGPVTD ALMDALDTPF TVFSPALPVN GRTVYQGYLF VMNQLLAESG MRHHPVNPMT DSYLPRLVEA QSTGRCGVVS AHVFEQGVDA VRQELARLQQ EGYRYAVLDA LTEHHLEIQG EALRDAPLVT GGSGLAIGLA RQWAQENGNQ ARKAGRPLAG RGVVLGSCS QMTNRQVAHY RQIAPAREVD VARCLSIETL AAYAHELAEW VLGQESVLAP LVFATASTDA LAAIQQQYGA QKASQAVETL FSQLAARLAA EGVTRFIVAG GETSGVVTQS LGIKGFHIGP TISPACRG
Source	Baculovirus
Target Names	otnK
Protein Names	Recommended name: Uncharacterized protein ygbK
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.