



# Recombinant Human Carboxypeptidase N subunit 2 (CPN2)

<b>Product Code</b>	CSB-BP005899HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P22792
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	CPMGCDCFV QEVFCSDEEL ATVPLDIPPY TKNIIFVETS FTTLETRAFG SNPNLTKVVF LNTQLCQFRP DAFGGLPRLE DLEVTGSSFL NLSTNIFSNL TSLGKLT LNF NMLEALPEGL FQHAALESL HLQGNQLQAL PRRLFQPLTH LKTLNLAQNL LAQLPEELFH PLTSLQTLKL SNNALSGLPQ GVFGKLGSLQ ELFLDSNNIS ELPPQVFSQL FCLERLWLQR NAITHLPLSI FASLGNLTF SLQWNMLRVL PAGLFAHTPC LVGLSLTHNQ LETVAEGTFA HLSNLRSLML SYNAITHLPA GIFRDLEELV KLYLGSNNLT ALHPALFQNL SKLELLSLSK NQLTTLPEGI FDTNYNLFNL ALHGPNWQCD CHLAYLFNWL QQYTDRLN QTYCAGPAYL KGQVVPALNE KQLVCPVTRD HLG FQVTWPD ESKAGGSWDL AVQERAARSQ CTYSNPEGTV VLACDQAQCR WLVNQLSPQQ GSLGLQYNAS QEWDLRSSCG SLRLTVSIEA RAAGP
<b>Source</b>	Baculovirus
<b>Target Names</b>	CPN2
<b>Protein Names</b>	Recommended name: Carboxypeptidase N subunit 2 Alternative name(s): Carboxypeptidase N 83 kDa chain Carboxypeptidase N large subunit Carboxypeptidase N polypeptide 2 Carboxypeptidase N regulatory subunit
<b>Expression Region</b>	22-545
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.