



Recombinant Zea mays Late embryogenesis abundant protein, group 3 (MGL3)

Product Code	CSB-MP673712ZAX
Abbreviation	MGL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q42376
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MASHQDKASY QAGETKARTE EKTGQAVGAT KDTAQHAKDR AADAAGHAAG KGQDAKEATK QKASDTGSYL GKKTDEAKHK AGETTEATKH KAGETTEAAK QKAGETTEAA KQKAGETTET TKQKAGETTE AARQKAADAM EAAKQKAAEA GQYAKDTAVS GKDKSGGVIQ QATEQVKSAA AGRKDAVMST LGMGGDNKQG DANTNTNTNT TKDSSTITRD H
Source	Mammalian cell
Target Names	MGL3
Protein Names	Recommended name: Late embryogenesis abundant protein, group 3 Short name= LEA
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.