



Recombinant Mono-ADP-ribosyltransferase C3 (c3)

Product Code	CSB-EP673832CEAI-B
Abbreviation	c3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46134
Product Type	Recombinant Protein
Immunogen Species	Hathewayia limosa (Clostridium limosum)
Purity	>85% (SDS-PAGE)
Sequence	PYADS FKEFTNIDEA RAWGDKQFAK YKLSSEKNA LTIYTRNAAR INGPLRANQG NTNGLPADIR KEVEQIDKSF TKMQTPENII LFRGDDPGYL GPDFENTILN RDGTINKAVF EQVKLRFK GK DRKEYGYIST SLVNGSAFAG RPIITKFKVL DGSKAGYIEP ISTFKGQLEV LLPRSSTYTI SDMQUIAPNNK QIITALLKR
Source	E.coli
Target Names	c3
Protein Names	Recommended name: Mono-ADP-ribosyltransferase C3 EC= 2.4.2.- Alternative name(s): Exoenzyme C3
Expression Region	46-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.