



Recombinant Escherichia coli Probable acetyl-CoA acetyltransferase (yqeF)

Product Code	CSB-EP673868ENV
Abbreviation	yqeF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46939
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKDVVIVGAL RTPIGCFRGA LAGHSAVELG SLVVKALIER TGVPAYAVDE VILGQVLTAG AGQNPARQSA IKGGLPNSVS AITINDVCGS GLKALHLATQ AIQCGEADIV IAGGQENMSR APHVLTDSTRT GAQLGNSQLV DSLVHDGLWD AFNDYHIGVT AENLAREYGI SRQLQDAYAL SSQQKARAAI DAGRFKDEIV PVMTQSNQQT LVVDTDEQPR TDASAEGLAR LNPSFDSLGS VTAGNASSIN DGAAAVMMMS EAKARALNLP VLARIRAFAS VGVDPALMGI APVYATRRCL ERVGWQLAEV DLIEANEFAA AQALSVGKML EWDERRVNVN GGAIALGHPI GASGCRILVS LVHEMVKRNA RKGLATLCIG GGQGVALTIE RDE
Source	E.coli
Target Names	yqeF
Protein Names	Recommended name: Probable acetyl-CoA acetyltransferase EC= 2.3.1.9 Alternative name(s): Acetoacetyl-CoA thiolase
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.