



# Recombinant Human Protein FAM188B (FAM188B)

<b>Product Code</b>	CSB-BP674576HU
<b>Abbreviation</b>	FAM188B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4G0A6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDSLFVEEVA ASLVREFLSR KGLKKTCTVM DQERPRSDLS INNRLDLRKY LHLEFLYKEN KAKENPLKTS LELITRYFLD HFGNTANNFT QDTPIPALSV PKKNNKVPSR CSETTLVNIY DLSDEDAGWR TSLSETSKAR HDNLDGDVLG NFVSSKRPPH KSKPMQTVPG ETPVLTSWE KIDKLHSEPS LDVKRMGENS RPKSGLIVRG MMSGPIASSP QDSFHRHYLR RSSPSSSSTQ PQEESRKVPE LFVCTQQDIL ASSNSSPSRT SLGQSELTV ERQKTTASSP PHLPSKRLPP WDRARPRDPS EDTPAVDGST DTDRLPLKLY LPGGNSRMTQ ERLERAFKRQ GSQPAPVRKN QLLPSDKVDG ELGALRLEDV EDELIREEVI LSPVPSVLKL QTASKPIDLS VAKEIKTLF GSSFCCFNEE WKLQSFSSFN TASLKYGIVQ NKGGPCGVLA AVQGCVLQKL LFEGDSKADC AQGLQPSDAH RTRCLVLALA DIVWRAGGRE RAVVALASRT QQFSPTGKYK ADGVLETTLT HSLTCYEDLV TFLQQSIHQF EVGPYGCILL TLSAILSRST ELIRQDFDVP TSHLIGAHGY CTQELVNLLL TGKAVSNVFN DVVELDSGDG NITLLRGIAA RSDIGFLSLF EHYNMCQVGC FLKTPRFPIW VVCSESHFSI LFSLQPGLLR DWRTERLFDL YYYDGLANQQ EQIRLTIDTT QTISEDNDND LVPPELCLIR TKWKGASVNW NGSDPIL
<b>Source</b>	Baculovirus
<b>Target Names</b>	MINDY4
<b>Protein Names</b>	Recommended name: Protein FAM188B
<b>Expression Region</b>	1-757
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.