



Recombinant Human Protein FAM188B (FAM188B)

Product Code	CSB-MP674576HU
Abbreviation	FAM188B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4G0A6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDSLFVEEVA ASLVREFLSR KGLKKTCTVM DQERPRSDLS INNRLDLRKY LHLEFLYKEN KAKENPLKTS LELITRYFLD HFGNTANNFT QDTPIPALSV PKKNNKVPSR CSETTLVNIY DLSDEDAGWR TSLSETSKAR HDNLDGDVLG NFVSSKRPPH KSKPMQTVPG ETPVLTSawe KIDKLHSEPS LDVKRMGENS RPKSGLIVRG MMSGPIASSP QDSFHRHYLR RSSPSSSSTQ PQEESRKVPE LFVCTQQDIL ASSNSSPSRT SLGQELSELV ERQKTTASSP PHLPSKRLPP WDRARPRDPS EDTPAVDGST DTDRLMPLKLY LPGGNSRMTQ ERLERAFKRQ GSQPAPVRKN QLLPSDKVDG ELGALRLEDV EDELIREEVI LSPVPSVLKL QTASKPIDLS VAKEIKTLF GSSFCCFNEE WKLQSF SFSN TASLKYGIVQ NKGGPCGVLA AVQGCVLQKL LFEGDskADC AqGLQPSDAH RTRCLVLALA DIVWRAGGRE RAVVALASRT QQFSPTGKYK ADGVLETTLT HSLTCYEDLV TFLQQSIHQF EVGPYGCILL TLSAILSRST ELIRQDFDVP TSHLIGAHGY CTQELVNLLL TGKAVSNVFN DVVELDsgDG NITLLRGIAA RSDIGFLSLF EHYNMCQVGC FLKTPRFPIW VVCSESHFSI LFSLQPGLLR DWRTERLFDL YYYDGLANQQ EQIRLTIDTT QTISEDtdND LVPPLELCIR TKWKGASVNW NGSDPIL
Source	Mammalian cell
Target Names	MINDY4
Protein Names	Recommended name: Protein FAM188B
Expression Region	1-757
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.