



Recombinant Zea mays Homeobox protein rough sheath 1 (RS1)

Product Code	CSB-BP674819ZAX
Abbreviation	RS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q41853
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MDQSFGNLGA GAGSSSSGGSN SKAAATAVSS SSFLQLPLST ASPAYYGAPL ALLHHAAAAP SSSQQHQQQQ HHHHYARHGA EMSAAEAEAI KAKIVAH PQY SALLAAYLDC QKVGAPPDVL ERLTAMAAKL DASAAGRHEP RDPELDQFME AYCNMLVKYR EELTRPIDEA MEFLKRVEAQ LDCISGGGGS SSARLSLADG KSEGVGSSSED DMDPNGREND PPEIDPRAED KELKYQLLKK YSGYLSSLRQ EFSKKKKKGK LPKEARQKLL HWWELHYKWP YPSETEKIAL AESTGLDQKQ INNWFINQRK RHWKPSDMP FVMMEGFHPQ NAAALYMDGP FMRDGM YRLG S
Source	Baculovirus
Target Names	RS1
Protein Names	Recommended name: Homeobox protein rough sheath 1
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.