



Recombinant Arabidopsis thaliana Glutamine synthetase, chloroplastic/mitochondrial (GLN2)

Product Code	CSB-YP674854DOA
Abbreviation	GLN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q43127
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	GFRVL ALQSDNSTVN RVETLLNLDLTPKPYSDRIIAE YIWIGGSGID LRSKSRTIEK PVEDPSELPK WNYDGSSTGQ APGEDSEVIL YPQAIFRDPF RGGNNILVIC DTWTPAGEPI PTNKRAKAAE IFSNKKVSGE VPWFGIEQEY TLLQQNVKWP LGWPVGAFPG PQGPYYCGVG ADKIWGRDIS DAHYKACLYA GINISGTNGE VMPGQWEFQV GPSVGIDAGD HVWCARYLLE RITEQAGVVL TLDPKPIEGD WNGAGCHTNY STKSMREEGG FEVIKKAILN LSLRHKEHIS AYGEGNERRL TGKHETASID QFSWGVANRG CSIRVGRDTE AKGKGYLEDR RPASNMDPYI VTSLLAETTL LWEPTLEAEA LAAQKLSLNV
Source	Yeast
Target Names	GLN2
Protein Names	Recommended name: Glutamine synthetase, chloroplastic/mitochondrial EC=6.3.1.2 Alternative name(s): GS2 Glutamate--ammonia ligase
Expression Region	46-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.