



# Recombinant Arabidopsis thaliana Peroxidase 64 (PER64)

<b>Product Code</b>	CSB-EP674875DOA-B
<b>Abbreviation</b>	PER64
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q43872
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LSPHYDDH TCPQADHIVT NAVKKAMSND QTVPAALLRM HFHDCFVRGC DGSVLLDSKG KNKAEKDGPP NISLHAFYVI DNAKKALEEQ CPGIVSCADI LSLAARDAVA LSGGPTWAVP KGRKDGRISK AIETRQLPAP TFNISQLRQN FGQRGLSMHD LVALSGGHTL GFAHCSSFQN RLHKFNTQKE VDPTLNPSFA ARLEGVCPAH NTVKNAGSNM DGTVTSFDNI YYKMLIQGKS LFSSDESLLA VPSTKKLVAK YANSNEEFER AFVKSMIKMS SISGNGNEVR LNCRRVR
<b>Source</b>	E.coli
<b>Target Names</b>	PER64
<b>Protein Names</b>	Recommended name: Peroxidase 64 Short name= Aterox P64 EC= 1.11.1.7 Alternative name(s): ATP17a PRXR4
<b>Expression Region</b>	23-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.