



# Recombinant *Drosophila melanogaster* Protein yorkie (yki)

<b>Product Code</b>	CSB-MP674941DLU
<b>Abbreviation</b>	yki
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q45VV3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCACLIAKII LCSFRLYTIS AFYMLTTMSA SSNTNSLIEK EIDDEDMLSP IKSNLNVVRV NQDTDDNLQA LFDSVLNPGD AKRPLQLPLR MRKLPNSFFT PPAPSHSRAN SADSTYDAGS QSSINIGNKA SIVQQPDGQS PIAAIPQLQI QPSPQHSRLA IHHSRARSSP ASLQQNYNVR ARSDAAAANN PNANPSSQQQ PAGPTFPENS AQEFPSGAPA SSAIDLAMN TCMSQDIPMS MQTVHKKQRS YDVISPIQLN RQLGALPPGW EQAKTNDGQI YYLNHTTKST QWEDPRIQYR QQQQLMAER IKQNDVLQTT KQTTTSTIAN NLGPLPDGWE QAVTESGDLY FINHIDRTTS WNDPRMQSGL SVLDCPDNLV SSLQIEDNLC SNLFNDAQAI VNPPSSHKPD DLEWYKIN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	yki
<b>Protein Names</b>	Recommended name: Protein yorkie
<b>Expression Region</b>	1-418
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.