



# Recombinant Escherichia coli Uncharacterized protein ygfT (ygfT)

<b>Product Code</b>	CSB-BP674973ENV
<b>Abbreviation</b>	ygfT
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q46820
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNKFIAAEAA ECIGCHACEI ACAVAHNQEN WPLSHSDFRP RIHVVGKGQA ANPVACHHCN NAPCVTACPV NALTFQSDSV QLDEQKCIGC KRCAIACPFV VVEMVDTIAQ KCDLCNQRSS GTQACIEVCP TQALRLMDDK GLQQIKVARQ RKTAAGKASS DAQPSRSAAL LPVNSRKGAD KISASERKTH FGEIYCGLDP QQATYESDRC VYCAEKANCN WHCPLHNAIP DYIRLVQEGK IIEAAELCHQ TSSLPEICGR VCPQDRLCEG ACTLKDHS GA VSIGNLERYI TDTALAMGWR PDVSKVVPRS EKVAVIGAGP AGLGCADILA RAGVQVDVFD RHPEIGGMLT FGIPPFKLDK TVLSQRREIF TAMGIDFHLN CEIGRDITFS DLTSEYDAVF IGVGTYGMMR ADLPHEAPG VIQALPFLTA HTRQLMGLPE SEEYPLTDVE GKRVVVLGGG DTTMDCLRTS IRLNAASVTC AYRRDEV SMP GSRKEVVNAR EEGVEFQFNV QPQYIACDED GRLTAVGLIR TAMGEPGPDG RRRPRPVAGS EFELPADVLI MAFGFQAHAM PWLQGSGIKL DKWGLIQTGD VGYLPTQTHL KKV FAGGDAV HGADLVVTAM AAGRQAARDM LTLFDTKAS
<b>Source</b>	Baculovirus
<b>Target Names</b>	ygfT
<b>Protein Names</b>	Recommended name: Uncharacterized protein ygfT
<b>Expression Region</b>	1-639
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .