




Recombinant Sheep Tumor necrosis factor (TNF), partial, partial

Product Code	CSB-YP023955SH
Relevance	Cytokine that binds to TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. It is mainly secreted by macrophages and can induce cell death of certain tumor cell lines. It is potent pyrogen causing fever by direct action or by stimulation of interleukin-1 secretion and is implicated in the induction of cachexia, Under certain conditions it can stimulate cell proliferation and induce cell differentiation. The TNF intracellular domain (ICD) form induces IL12 production in dendritic cells.
Abbreviation	Recombinant Sheep TNF protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23383
Alias	Cachectin;TNF-alphaTumor necrosis factor ligand superfamily member 2 ;TNF-a
Product Type	Recombinant Protein
Immunogen Species	Ovis aries (Sheep)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	LRSSSQASNNKPVAHVVANISAPGQLRWGDSYANALMANGVELKDNQLVVPT DGLYLIYSQVLFRGHGCPSTPLFLTHTISRIAVSYQTKVNILSAIKSPCHRETLEG AEAKPWYEPIYQGGVFQLEKGDRLSAEINLPEYLDYAESGQVYFGIIAL
Source	Yeast
Target Names	TNF
Protein Names	Cachectin;TNF-alphaTumor necrosis factor ligand superfamily member 2 ;TNF-a
Expression Region	78-234aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	19.2 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.