



Recombinant Escherichia coli Copper resistance protein A (pcoA)

Product Code	CSB-BP675105ENL
Abbreviation	pcoA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q47452
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	LQGTQFDLT IGETAVNITG SERQAKTING GLPGPVLRWK EGDTITLKVK NRLNEQTSIH WHGIILPANM DGVPGLSFMG IEPDDTYVYT FKVKQNGTYW YHSHSGLQEQ EGVYGAIID AREPEPFAYD REHVVMLSDW TDENPHSLLK KLKKQSDYYN FNKPTVGSFF RDVNTRGLSA TIADRKMWAE MKMNPTDLAD VSGYTYTYLM NGQAPLKNWT GLFRPGEKIR LRFINGSAMT YFDIRIPGLK MTVVAADGQY VNPVTVDEFR IAVAETYDVI VEPQGEAYTI FAQSMDRGTGY ARGTLATREG LSAAVPPLDP RPLLTMEDMG MGGMGHDMAG MDHSQMGGMD NSGEMMSMDG ADLPDSGTSS APMDHSSMAG MDHSRMAGMP GMQSHPASET DNPLVDMQAM SVSPKLNDPG IGLRNNGRKV LTYADLKSFR EDPDGREPGR TIELHLTGHM EKFAWSFNIGI KFSDAAPVLL KYGERLRITL INDTMMTHPI HLHGMWSDLE DENGFMVRK HTIDVPPGTK RSYRVTADAL GRWAYHCHLL YHMEMGMFRE VRVEE
Source	Baculovirus
Target Names	pcoA
Protein Names	Recommended name: Copper resistance protein A
Expression Region	42-605
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.