



Recombinant Escherichia coli Taurine-binding periplasmic protein (tauA)

| | |
|--------------------------|---|
| Product Code | CSB-YP675110ENV |
| Abbreviation | tauA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q47537 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | VNVTVAYQ TSAEPAKVAQ ADNTFAKESG ATVDWRKFDS GASIVRALAS GDVQIGNLGS SPLAVAASQQ VPIEVFLLAS KLGNSEALVV KKTISKPEDL IGKRIAVPFI STTHYSLLA LKHGWIKPGQ VEIVNLQPPA IIAAWQRGDI DGAYVWAPAV NALEKDGKVL TDSEQVGQWG APTLDVWVVR KDFAEKHPEV VKAFAKSAID AQQPYIANPD VWLKQPENIS KLARLSGVPE GDVPGLVKGN TYLTPQQQTA ELTGPNKAI IDTAQFLKEQ GKVPVANDY SQYVTSRFVQ |
| Source | Yeast |
| Target Names | tauA |
| Protein Names | Recommended name: Taurine-binding periplasmic protein Alternative name(s): Sulfate starvation-induced protein 1 Short name= SSI1 |
| Expression Region | 23-320 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |