



Recombinant Human Myogenic factor 6 (MYF6)

Product Code	CSB-MP015288HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P23409
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMMDLFETGS YFFYLDGENV TLQPLEVAEG SPLYPGSDGT LSPCQDQMPP EAGSDSSGEE HVLAPPGLQP PHCPGQCLIW ACKTCKRKSA PTDRRKAATL RERRRLKKIN EAFEALKRRT VANPNQRLPK VEILRSAISY IERLQDLLHR LDQQEKMQEL GVDPFYSYRPK QENLEGADFL RTCSSQWPSV SDHSRGLVIT AKEGGASIDS SASSSLRCLS SIVDSISSEE RKLPCVEEVV EK
Source	Mammalian cell
Target Names	MYF6
Protein Names	Recommended name: Myogenic factor 6 Short name= Myf-6 Alternative name(s): Class C basic helix-loop-helix protein 4 Short name= bHLHc4 Muscle- specific regulatory factor 4
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.