



Recombinant Mouse Protein FAM5C (Fam5c)

Product Code	CSB-MP675462MO
Abbreviation	Fam5c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q499E0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AVSDQHA TSPFDWLLSD KGPFHRSQEY TDFVDRSRQG FSTRYKIYRE FGRWKVNNLA VERRNFLGSP LPLAPEFFRN IRLLGRRPTL QQITENLIKK YGTHFLLSAT LGGEESLTIF VDKRKLSKRS EGSETSTNSS SVTLETLHLQ AASYFIDRDS TLRLHHIQI ASTAIKVET RTGPLGCSNY DNLDVSSVL VQSPENKIQL QGLQVLLPDY LQERFVQAAL SYIACNSEGE FICKDNDWC HCGPKFPECN CPSMDIQAME ENLLRITETW KAYNSDFEDS DEFKFFMKRL PMNYFLNTST IMHLWTMDSN FQRRYELEN SMKQLFLKAH RIVHKLFSLS KRCHKQPLIS LPRQRTSTYW LTRIQSFLYC NENGLLGSFS EETHSCTCPN DQVVCTAFLP CTVGDAACL TCAPDNRTRC GTCNTGYMLS QGLCKPEVAE STDHYIGFET DLQDLEMKYL LQKTDRRIEV HAIFISNDRM LNSWFDPSWR KRMLLTLKSN KYKSSLVHMI LGLSLQICLT KNSTLEPVLA VYVNPFGGSH SESWFMPVSE SSFPDWERTK LDLPLQCYNW TLTLGNKWKT FFETVHIYLR SRIKANGPNS NESIYYEPL FIDPSRNLGY MKINNIQVFG YSMHFDPEAI RDLILQLDYP YTQGSQDSAL LQLEIRDRV NKLSPPGQRR LDLFSCLLRH RLKLSTSEVV RIQSALQAFN AKLPNTVDYD TTKLCS
Source	Mammalian cell
Target Names	Brinp3
Protein Names	Recommended name: Protein FAM5C
Expression Region	34-766
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.