



# Recombinant *Oryza sativa* subsp. *japonica* Tubulin beta-3 chain (TUBB3)

<b>Product Code</b>	CSB-YP675920OFG
<b>Abbreviation</b>	TUBB3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q40665
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MREILHIQGG QCGNQIGAKF WEVICDEHGI DHTGKYSGDS DLQLERINVY YNEASGGRYV PRAVLMDLEP GTMDSVRS GP YGQIFRPDNF VFGQSGAGNN WAKGHYTEGA ELIDSVLDVV RKEAENCDCDCL QGFQVCHSLG GGTGSGMGTL LISKIREEYP DRMMLTFSVF PSPKVS DTVV EPYNATLSVH QLVENADECM VLDNEALYDI CFRTLKLATP TFGDLNHLIS ATMSGVTCCL RFPGQLNSDL RKLAVNLIPF PRLHFFMVGF APLTSRGSQQ YRALTVPELT QQMWDANKMM CAADPRHGRY LTASAMFRGK MSTKEVDEQM LNVQNKSSY FVEWIPNNVK SSSVCDIPPIG LKMASTFIGN STSIQEMFRR VSEQFTAMFR RKAFLHWYTG EGMDEMEFTE AESNMNDLVA EYQQYQDATA EEEDYEEEEEE DEEVAA
<b>Source</b>	Yeast
<b>Target Names</b>	TUBB3
<b>Protein Names</b>	Recommended name: Tubulin beta-3 chain Alternative name(s): Beta-3-tubulin
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.