



Recombinant Solanum tuberosum 14-3-3-like protein

Product Code	CSB-BP675933FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q41418
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MADSREENVY MAKLAEQAER YEEMVEFMEK VAKVDVEELT VEERNLLSVA YKNVIGARRA SWRIISSIEQ KEESRGNEHD VSSIKEYRVK IEAELSKICD GILSLLESHL VPSASTAESK VFYLKMKGDY HRYLAEFKGT ARERKEAAEN TLLAYKSAQD IALAELTPTH PIRLGLALNF SVFYEILNS PDRACNLAKQ AFDDAIAELD TLGEESYKDS TLIMQLLRDN LTLWTSDDTD DAEDEIREGS KQESGDGQQ
Source	Baculovirus
Protein Names	Recommended name: 14-3-3-like protein
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.