



# Recombinant Mouse Plasma serine protease inhibitor (Serpina5)

<b>Product Code</b>	CSB-EP021061MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P70458
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SKKKKA KESSVGAVGP PSSKDFAFRL YRALVSESPG QNVFFSPLSV SMSLGMLSLG AGLKTKTQIL DGLGLSLQQG QEDKLHKGfQ QLLQRFRQPS DGLQLSLGSA LFKDPAVHIR DDFLSAMKTL YMSDTFSTNF GNPEIAKKQI NNYVAKQTKG KIVDFIKDLD STHVMIVVNY IFFKAKWQTA FSETNTHKMD FHVTPKRTTQ VPMMNREDGY SYLDQNIsc TVVGIPYQGN AIALFILPSE GKMKQVEDGL DERTLRNWLK MFTKRRLDLY LPKFSIEATY KLENVLPKLG IQDVFTTHAD LSGITDHTNI KLSEMVHKSM MEVEESGTTA AAITGAIFTF RSARPSSLKI EFTRPFLTL MEDSHILFVG KVTRP
<b>Source</b>	E.coli
<b>Target Names</b>	Serpina5
<b>Protein Names</b>	Recommended name: Plasma serine protease inhibitor Alternative name(s): Plasminogen activator inhibitor 3 Short name= PAI-3 Short name= PAI3 Protein C inhibitor Short name= PCI Serpin A5
<b>Expression Region</b>	25-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.