



Recombinant Bradyrhizobium japonicum Probable farnesyl diphosphate synthase (fppS)

Product Code	CSB-EP676034BVW-B
Abbreviation	fppS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q45220
Product Type	Recombinant Protein
Immunogen Species	Bradyrhizobium diazoefficiens (strain JCM 10833 / IAM 13628 / NBRC 14792 / USDA 110)
Purity	>85% (SDS-PAGE)
Sequence	MQTGSTPHDD RAGVSANGIL GAQARGASGR LLPEIWMQDG AKRVEEALAR LLCAEDDGET ELMAAMRYAT LHGGKRTRAL LCLAAGALAD TPAHMLDDVG AAIEMMHACT LVHDDLPA MD DDVLRRLPT VHVKFG EATA ILVGDALQAH AFLTLASLNA PGDSPIALVR ELAQAVSAEG AAGGQAIDLS LVGKHVELDR IVAMHRMKSG ALVRASVRMG ALCAVGVNAA HAALYCALDH YSACFGLALQ VIDDILDVTA DTAALGKTPG KDAAAQKPTC ASIMGLQEAR QFALDLLRDA GEAIAPLGPR AERLAQLIQR ANAYLFKHAP RA
Source	E.coli
Target Names	fppS
Protein Names	Recommended name: Probable farnesyl diphosphate synthase Short name= FPP synthase EC= 2.5.1.10 Alternative name(s): (2E,6E)-farnesyl diphosphate synthase Geranyltranstransferase
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.