



Recombinant Human Poly [ADP-ribose] polymerase 15 (PARP15)

Product Code	CSB-EP676066HU-B
Abbreviation	PARP15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q460N3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MHGVAGVTSR AGRDREAGSM LPAGNRGARK ASRRSSSRSM SRDNKFSKKD CLSIRNVVAS IQTKEGLNLK LISGDVLYIW ADVIVNSVPM NLQLGGGPLS RAFLQKAGPM LQKELDDRRR ETEEKVGNIF MTSGCNLDCK AVLHAVAPYW NNGAETSWQI MANIIKKCLT TVEVLSFSSI TFPMIGTGSL QFPKAVFAKL ILSEVFEYSS STRPITSPLQ EVHFLVYTND DEGCQAFLE FTNWSRINPN KARIPMAGDT QGVVGTVSKP CFTAYEMKIG AITFQVATGD IATEQVDVIV NSTARTFNKR SGVSRAILEG AGQAVESECA VLAAQPHRDF IITPGGCLKC KIIHVPGGK DVRKTVTSVL EECEQRKYTS VSLPAIGTGN AGKNPITVAD NIIDAIVDFS SQHSTPSLKT VKVVIFQPEL LNIFYDSMKK RDLASASLNFQ STFSMTTCNL PEHWTDNMNHQ LFCMVQLEPG QSEYNTIKDK FTRTCSSYAI EKIERIQNAF LWQSYQVKKR QMDIKNDHKN NERLLFHGTD ADSVPYVNQH GFNRSCAGKN AVSYGKGYF AVDASYSKD TYSKPDSNGR KHMVVRVLT GVFTKGRAGL VTPPPKNPHN PTDLFDSVTN NTRSPKLFVV FFDNQAYPEY LITFTA
Source	E.coli
Target Names	PARP15
Protein Names	Recommended name: Poly [ADP-ribose] polymerase 15 Short name= PARP-15 EC= 2.4.2.30 Alternative name(s): B-aggressive lymphoma protein 3
Expression Region	1-656
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.