



# Recombinant Human Poly [ADP-ribose] polymerase 15 (PARP15)

<b>Product Code</b>	CSB-MP676066HU
<b>Abbreviation</b>	PARP15
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q460N3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MHGVAGVTSR AGRDREAGSM LPAGNRGARK ASRRSSSRSM SRDNKFSKKD CLSIRNVVAS IQTKEGLNLK LISGDVLYIW ADVIVNSVPM NLQLGGGPLS RAFLQKAGPM LQKELDDRRR ETEEKVGNIF MTSGCNLDCK AVLHAVAPYW NNGAETSWQI MANIIKKCLT TVEVLSFSSI TFPMIGTGSL QFPKAVFAKL ILSEVFEYSS STRPITSPLQ EVHFLVYTND DEGCQAFLE FTNWSRINPN KARIPMAGDT QGVVGTVSKP CFTAYEMKIG AITFQVATGD IATEQVDVIV NSTARTFNK SGVSRAILEG AGQAVESECA VLAAQPHRDF IITPGGCLKC KIIHVPGGK DVRKTVTSVL EECEQRKYTS VSLPAIGTGN AGKNPITVAD NIIDAIVDFS SQHSTPSLKT VKVVIFQPEL LNIFYDSMKK RDLASLNFQ STFSMTTCNL PEHWTDNMNHQ LFCMVQLEPG QSEYNTIKDK FTRTCSSYAI EKIERIQNAF LWQSYQVKKR QMDIKNDHKN NERLLFHGTD ADSVPYVNQH GFNRSCAGKN AVSYGKGYF AVDASYSKAD TYSKPDSNGR KHMVVRVLT GVFTKGRAGL VTPPPKNPHN PTDLFDSVTN NTRSPKLFVV FFDNQAYPEY LITFTA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PARP15
<b>Protein Names</b>	Recommended name: Poly [ADP-ribose] polymerase 15 Short name= PARP-15 EC= 2.4.2.30 Alternative name(s): B-aggressive lymphoma protein 3
<b>Expression Region</b>	1-656
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.