



Recombinant Escherichia coli 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase (ispD)

Product Code	CSB-MP676100ENV
Abbreviation	ispD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46893
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	ATTHLDVCA VVPAAGFGRR MQTECPKQYL SIGNQTILEH SVHALLAHPR VKRVVIAISP GDSRFAQLPL ANHPQITVVD GGDERADSVL AGLKAAGDAQ WVLVHDAARP CLHQDDLARL LALSETSRTG GILAAPVRDT MKRAEPGKNA IAHTVDRNGL WHALTPQFFP RELLHDCLTR ALNEGATITD EASALEYCGF HPQLVEGRAD NIKVTRPEDL ALAEFYLTRT IHQENT
Source	Mammalian cell
Target Names	ispD
Protein Names	Recommended name: 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase EC= 2.7.7.60 Alternative name(s): 4-diphosphocytidyl-2C-methyl-D-erythritol synthase CDP-ME synthase MEP cytidyltransferase Short name= MCT
Expression Region	2-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.