



# Recombinant Colwellia psychrerythraea DNA polymerase IV (dinB)

<b>Product Code</b>	CSB-MP676401CBAB
<b>Abbreviation</b>	dinB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q487H6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Colwellia psychrerythraea (strain 34H / ATCC BAA-681) (Vibrio psychroerythus)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGNQKKIHI DMDCFYAAIE MRDFPEYQNI PLAVGGDGPR SVLCTSNYQA RQFGVRSAMP AIKAKQLCPH LKIVHGRMDV YKETSKNIRE IFSRYTDLIE PLSLDEAYLD VTDATMCQGS ATLIAERIRA DIFNELNLTA SAGIAPNKFL AKIASDENKP NGQCVITPDK VANFVEQLSL KKIPIGIGPKT FEKLNHRHGYV TCADVRQSNI RALQNIVGKF ANSLYLKSHG VDNRDLEVSR QRKSLAIETT LAHDISTQDE CKLVIDSLYQ KLLTRLAPHS NREIIRQGVK LKFTDFNQTT VETQSNECQQ ALFISLLSKA YSRSNKRGVR LVGLTLGFAD SPGESQQLSL SL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	dinB
<b>Protein Names</b>	Recommended name: DNA polymerase IV Short name= Pol IV EC= 2.7.7.7
<b>Expression Region</b>	1-352
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.