



Recombinant Rat DnaJ homolog subfamily C member 10 (Dnajc10)

Product Code	CSB-EP676589RA-B
Abbreviation	Dnajc10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q498R3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	DQNFYSLL GVSKTASSRE IRQAFKKLAL KLHPDKNPNN PNAHGDFLKI NRAYEVLKDE DLRKKYDKYG EKGLEDNQGG QYESWSYYRY DFGIYDDDPE IITLERREFD AAVNSGELWF VNFYSPGCSH CHDLAPTWRE FAKEVDGLLR IGAVNCGDDR MLCRMKGVNS YPSLFIFRSG MAAVKYNGDR SKESLVSFAM QHVRTTVTEL STGNFVNAIE TAAAGIGWL ITFCFKGEDC LTPQTRLRLS GMLDGLVNVG WVDCDTQDSL CKSLDATAST TAYFPPGATL NNKEKSSVLF LNSLDAKEIY MEIIHNLPDF ELLSANKLED RLAHHRWLVF FHFQKNNAN DPELKKLKTLLKNEHIQVGR FDCSSAPGIC SDLYVFQSCS AVFKGQGTKE YEIHHGKKIL YDILAFKES VNSHVTTLGP QNFPASDKEP WLVDFFAPWC PPCRALLPEL RCASTLLYGQ LKVGTLDCIT HEGLCNMYNI QAYPTTVVFN QSSVHEYEGH HSAEQILEFI EDLRNPSVVS LPTTFNELV KQRKHDEVWM VDFYSPWCHP CQVLMPEWKR MARTLTGLIN VGSVDCQQYH SFCTQENVQR YPEIRFYPQK SSRAYQYHSY NGWNRDAYSL RSWGLGFLPQ ASIDLTPQTF NEKVLQGKTH WVIDFYAPWC GPCQNFAPF ELLARMIK GK VKAGKVDCQA YPQTCQKAGI RAYPSVKLYL YERAKKSIWE EQINSRDAKT IAALIYGKLE TFQSQVKNK DEL
Source	E.coli
Target Names	Dnajc10
Protein Names	Recommended name: DnaJ homolog subfamily C member 10
Expression Region	33-793
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.