



# Recombinant Rat DnaJ homolog subfamily C member 10 (Dnajc10)

<b>Product Code</b>	CSB-YP676589RA
<b>Abbreviation</b>	Dnajc10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q498R3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	DQNFYSLL GVSKTASSRE IRQAFKKLAL KLHPDKNPNN PNAHGDFLKI NRAYEVLKDE DLRKKYDKYG EKGLEDNQGG QYESWSYYRY DFGIYDDDPE IITLERREFD AAVNSGELWF VNFYSPGCSH CHDLAPTWRE FAKEVDGLLR IGAVNCGDDR MLCRMKGVNS YPSLFIFRSG MAAVKYNGDR SKESLVSFAM QHVRTTVEL STGNFVNAIE TAAAGIGWL ITFCFKGEDC LTPQTRLRLS GMLDGLVNVG WVDCDTQDSL CKSLDATAST TAYFPPGATL NNKEKSSVLF LNSLDAKEIY MEIIHNLPDF ELLSANKLED RLAHHRWLVF FHFQKNNAN DPELKKLKTLLKNEHIQVGR FDCSSAPGIC SDLYVFQSCS AVFKGQGTKE YEIHHGKKIL YDILAFKES VNSHVTTLGP QNFPASDKEP WLVDFFAPWC PPCRALLPEL RCASTLLYGQ LKVGTLDCIT HEGLCNMYNI QAYPTTVVFN QSSVHEYEGH HSAEQILEFI EDLRNPSVVS LPTTFNELV KQRKHDEVWM VDFYSPWCHP CQVLMPEWKR MARTLTGLIN VGSVDCQQYH SFCTQENVQR YPEIRFYPQK SSRAYQYHSY NGWNRDAYSL RSWGLGFLPQ ASIDLTPQTF NEKVLQGKTH WVIDFYAPWC GPCQNFAPF ELLARMIK GK VKAGKVDCQA YPQTCQKAGI RAYPSVKLYL YERAKKSIWE EQINSRDAKT IAALIYGKLE TFQSQVKNK DEL
<b>Source</b>	Yeast
<b>Target Names</b>	Dnajc10
<b>Protein Names</b>	Recommended name: DnaJ homolog subfamily C member 10
<b>Expression Region</b>	33-793
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.