



# Recombinant Human Vacuolar protein sorting-associated protein 26B (VPS26B)

<b>Product Code</b>	CSB-EP676824HU-B
<b>Abbreviation</b>	VPS26B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4G0F5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFFGFGQSV EVEILLNDAE SRKRAEHKTE DGKKEKYFLF YDGETVSGKV SLALKNPNKR LEHQGIKIEF IGQIELYYDR GNHHEFVSLV KDLARPGEIT QSQAQDFEFT HVEKPYESYT GQNVKLRFLYFL RATISRRLND VVKEMDIVVH TLSTYPELNS SIKMEVGIED CLHIEFEYNK SKYHLKDVIV GKIYFLLVRI KIKHMEIDII KRETTGTGPN VYHENDTIAK YEIMDGAPVR GESIPIRLFL AGYELTPTMR DINKKFSVRY YLNLVLIDEE ERRYFKQQEV VLWRKGDIVR KSMHQAAIA SQRFEGTSL GEVRTPSQLS DNPCRQ
<b>Source</b>	E.coli
<b>Target Names</b>	VPS26B
<b>Protein Names</b>	Recommended name: Vacuolar protein sorting-associated protein 26B Alternative name(s): Vesicle protein sorting 26B
<b>Expression Region</b>	1-336
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.