



Recombinant Arabidopsis thaliana 60S ribosomal protein L8-3 (RPL8C)

Product Code	CSB-EP677072DOA
Abbreviation	RPL8C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q42064
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRVIRAQRK GAGSVFKSHT HHRKGPAPFR SLDFGERNGY LKGVVTEIIH DPGRGAPLAR VAFRHPFRFK KQKELFVAEE GMYTGQFLYC GKKATLVVGN VLPLRSIPEG AVICNVEHHV GDRGVFARAS GDYAIVIAHN PDNDTSRIKL PSGSKKIVPS GCRAMIGQVA GGGRTEKPML KAGNAYHKYR VKRNCWPKVR GVAMNPVEHP HGGGNHQHIG HASTVRRDAP PGKKVGLIAA RRTGRLRGQA AALASKAD
Source	E.coli
Target Names	RPL8C
Protein Names	Recommended name: 60S ribosomal protein L8-3
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.