



Recombinant *Gossypium hirsutum* Acyl-[acyl-carrier-protein] desaturase, chloroplastic

Product Code	CSB-BP677090GHB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q42770
Product Type	Recombinant Protein
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Purity	>85% (SDS-PAGE)
Sequence	ISTIPSG SKEVGNLKKP FTPPKEVPVQ ITHSMPPHKI EIFKSLEGWA ENNILTHLKP VEKCWQPADF LPDPNSDGFH EQVKELRERA KEIPDDYFVV LVGDMITEEA LSTYQTMLNT LDGTRDETGA SLTPWAIWTR AWTAEENRHG DLLNKYLYLS GRVDMRQIER TIQYLIGSGM DPHTENSPYR GFIYTSFQER ATFISHGNTG RLAKEYGDIN LAQICGSIAS DEKRHETAYT KIVEKLFEID PDETVLAFAD MMKKKIAMPA EFIYDGRDYN LFDHYSAVAQ RIGVYTAKDY VDIVEHLVDR WKVKELAGLS AEGRKAQDYL CSLPSRIRRL EERAQEKAQG STPVSPFSWI FDREVKL
Source	Baculovirus
Protein Names	Recommended name: Acyl-[acyl-carrier-protein] desaturase, chloroplastic EC=1.14.19.2 Alternative name(s): Stearoyl-ACP desaturase
Expression Region	34-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.