



Recombinant *Desulfarculus baarsii* Desulfoferrodoxin (dfx)

Product Code	CSB-MP677203DBAJ
Abbreviation	dfx
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46495
Product Type	Recombinant Protein
Immunogen Species	<i>Desulfarculus baarsii</i> (strain ATCC 33931 / DSM 2075 / VKM B-1802 / 2st14)
Purity	>85% (SDS-PAGE)
Sequence	PERLQVYKC EVCGNIVEVL NGGIGELVCC NQDMKLMSSEN TVDAAKEKHHV PVIEKIDGGY KVKVGAHAHP MEEKHYIQWI ELLADDKCYT QFLKPGQAPE AVFLIEAAKV VAREYCNIHG HWKAEN
Source	Mammalian cell
Target Names	dfx
Protein Names	Recommended name: Desulfoferrodoxin Short name= Dfx EC= 1.15.1.2 Alternative name(s): Superoxide reductase Short name= SOR
Expression Region	2-126
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.