



Recombinant Escherichia coli Putative uncharacterized protein ygeQ (ygeQ)

Product Code	CSB-EP677216ENV-B
Abbreviation	ygeQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46797
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MVIYAFNKRL MEYFMKGKSA LLLLLAGIFS CGTCQATGAE VTSESVFNIL NSTGAATDKS YLSLNPDKYP NYRLLIHS AK LQNEIKSHYT KDEIQGLLTL TENTRKLTLT EKPWGTFILA STFEDDKTAA ETHYDAVWLR DSLWGYMALV SDQGN SVA AK K VLLTLWDYM STPDQIKRMQ DVISNPKRLD GIPVQMNAVH IRFDSNSPVM ADVQEEGKPQ LWNHKQNDAL GLYLDLLIQA INTGTINAED WQKGDRLKSV ALLIAYLDKA NFYVMEDS
Source	E.coli
Target Names	ygeQ
Protein Names	Recommended name: Putative uncharacterized protein ygeQ
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.