



Recombinant Escherichia coli Putative type II secretion system C-type protein YghF (yghF)

Product Code	CSB-MP677220ENV
Abbreviation	yghF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46834
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	AEYTAV SQPVNKPSRV DAKTFDKNDV QLISQQNWFG KYQPVAAQVK QPEPVPVAET RLNVVLRGIA FGARPGAVIE EGGKQQVYLQ GETLGSHNAV IEEINRDHVM LRYQGKIERL SLAEEERSTV AVTNKKAUSD EAKQAVAEP VSVPEIPAA VRQALAKDPQ KIFNYIQLTP VRKEGIVGYA AKPGADRSLF DASGFKEGDI AIALNQQDFT DPRAMIALMR QLPSMDSIQL TVLRKGARHD ISIALRYRLF CPPYWKKS
Source	Mammalian cell
Target Names	yghF
Protein Names	Recommended name: Putative type II secretion system C-type protein YghF Alternative name(s): Putative general secretion pathway C-type protein YghF
Expression Region	25-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.