



Recombinant *Dechloromonas aromatica* Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-EP677400DAAH
Abbreviation	psd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q47C25
Product Type	Recombinant Protein
Immunogen Species	<i>Dechloromonas aromatica</i> (strain RCB)
Purity	>85% (SDS-PAGE)
Sequence	MSDRLAVLPQ YLIPKQALTV LAGKLASAKA GGLTTSVIRW FVRRYVNMT EAANPDIASY KSFNEFFTRP LKDGARPAAD ADFLCPVDGA ISQYGTIDRD QIFQAKGHSY STTALVGGDR KLAEQFENGs FATLYLSPRD YHRIHMPCDG KLTRMIYVPG ALFSVNPTTA RGVPLFARN ERVICVFESE FGSFVLTLVG ATIVGSMATV WHGTINPPRP GVIREWRYDE QNIRLKKGQE MGRFLLG
Source	<i>E. coli</i>
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.