



Recombinant Trypanosoma cruzi Flap endonuclease 1 (FEN1)

Product Code	CSB-MP677857TAAG
Abbreviation	FEN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4DKQ5
Product Type	Recombinant Protein
Immunogen Species	Trypanosoma cruzi (strain CL Brener)
Purity	>85% (SDS-PAGE)
Sequence	MGILGLSKLL YDRSPAIRE RELKNYFGRR IAIDASMTIY QFIAMKGFQ DGQGMELTNE AGEVTSHLNG LFARTLRMVD EGLRPIYVFD GKPPTLKASE LQERRQRAEE AQLFDTAKE EGNDELMEKM SKRTVRVSRE QLEEAKKLLQ LMGIPVVQAP SEAEAQCAEL VKKKKAWAVA TEDMDALTFG APVMLRHLTY SEAKKRPIAE FHLDEILGIT GLTMTQFIDL CILLGCDYVP KIPGIGPQKA WEGIKKHGDI ETLLQSLDAG RHSVPEGFHY EEARQFFLKP EVTPAEEIEI QFREPEEGL VKFLVEEKLF NKDRV LKGIQ RLRNALTRKT QGRLDQFFTI TRPVTKPNTC DAKAGVKRGH SAIALSGTLQ QKGSSGHKKV VKK
Source	Mammalian cell
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.