



Recombinant *Oryza sativa* subsp. japonica Myb-related protein MYBAS2 (MYBAS2)

Product Code	CSB-EP678303OFG-B
Abbreviation	MYBAS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4JL76
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVTVREEIRK GPWTEQEDLQ LVCTVRLFGE RRWDFIAKVS GLNRTGKSCR LRWVNYLHPG LKRGRMSPHE ERLILELHAR WGNRWSRIAR RLPGRTDNEI KNYWRTHMRK KAQERKSNMS PSSSSSSSLTY QSCHPETPSM IIGIEEQELH GGSGCITSIM KSTPVDMDGY PMDQIWMEIE APNVLPGPCF DEAKDSASNS LSGPLLPYPM WDYYCPETCL RMDDEIKVAP QFGYGKGVGP CY
Source	E.coli
Target Names	MYBAS2
Protein Names	Recommended name: Myb-related protein MYBAS2
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.