



Recombinant Danio rerio Fructose-bisphosphate aldolase C-A (aldoca)

Product Code	CSB-MP679662DIL
Abbreviation	aldoca
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4KMC8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	THQFPPLTT EQKKELHEIA LRIVSPGKGI LAADESIGSM GKRLNQIGVE NNEENRRLFR QVLFTADDRI DNCIGGVIFF HETLYQNSDD GVPFVKMIKD KGITIGIKVD KGVVPLPGTN GETATQGLDG LSERCAQYKK DGADFAKWRC VMKISSETPS NLCITENAKV LARYASICQQ HGIVPIVEPE ILPDGDHNLK RCQFVTERVL AAVYKAMFDH HVYLEGTLK PNMVTPGHGC PTKYSAEEVA MATVTALRRT VPPAVTGVTF LSGGQSEEEA SINLSAINNC RLVKPWALTF SFGRALQASA LKTWRGQREN EAAATEEFIK RAEINSLASQ GKYTVCGDSS GATGLSHYLS SYAY
Source	Mammalian cell
Target Names	aldoca
Protein Names	Recommended name: Fructose-bisphosphate aldolase C-A EC= 4.1.2.13 Alternative name(s): Aldolase C-like Brain-type aldolase-A
Expression Region	2-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.