



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-D4-2 (CYCD4-2)

Product Code	CSB-BP679765OFG
Abbreviation	CYCD4-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4KYM5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAPSSSSCHD AAASMLLCAE DNSSILWLED EEGEVGERRS GGCRSMVGDL AAGGGGGSGG GGVEEEEDMF PRQSEECVAS LVEREQAHMP RADYGERLRG GGGDVDLRVR SEAIGWIWEV YTYYNFSSVT AYLAVNYLDR FLSQYELPEG RDWMTQLLSV ACLSIAAKME ETVVPQCLDL QIGEPFLFE VETIHRMELL VLTNLNWRMQ AVTPFSYIDY FLRKLNSGNA APRSWLLRSS ELILRIAAGT GFLEFRPSEI AAAVAATVAG EATGVVEEDI AEAFTHVDKG RVLQCQEAIQ DHHYSMATIN TVQPKPASTR RGSASASSSS VPESPAVLD AGCLSYKSDD TDAATIASHG GGRRKSCFDS SPVTSKKRRK LSR
Source	Baculovirus
Target Names	CYCD4-2
Protein Names	Recommended name: Cyclin-D4-2 Alternative name(s): G1/S-specific cyclin-D4-2 Short name= CycD4;2
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.