



Recombinant Danio rerio tRNA selenocysteine 1-associated protein 1-like (trnau1apl)

Product Code	CSB-MP679871DIL
Abbreviation	trnau1apl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4KM14
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MFNRMTSLWM GDLDPYMDEN FIKQAFSTMG ETAFGVKIIT HRVTGGSAGY CFVEMADEAS VDRCVQLNG KLVPGSNPPR KFKLNYATYG KRPEPGPEFS VFGDLTSEV DDYQLHQFFL KKFPSCKGAK VVTDPYGNSR GYGFVKFSDE NEQKKALEEF QNASGLGGKP IRISIAVNKG NKASTYHNQN NTYNTNYQQQ YYRQPNSYY PQWGYDQYSG YNYGYNPYAA PPPMMGPPPP MGMPMPPPDM QGSTEAHDGT EEVEEDPSED PNPQVDVEEL NRQYMERSEE LYDSLMECHW LPMDTITSDI SMMNGS
Source	Mammalian cell
Target Names	trnau1apl
Protein Names	Recommended name: tRNA selenocysteine 1-associated protein 1-like Alternative name(s): tRNA selenocysteine 1-associated protein 1 tRNA selenocysteine-associated protein 1
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.