



Recombinant Mouse Tyrosine-protein kinase Tec (Tec)

Product Code	CSB-EP023367MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24604
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p>MNFNTILEEI LIKRSQQKKK TSELLNYKERL CVLPKSVLSY YEGRAEKKYR KGVIDISKIK CVEIVKNDDG VIPCQNKFPF QVVHDANTLY IFAPSPQSRD RWVKKLKEEI KNNNNIMIKY HPKFWADGSY QCCRQTEKLA PGCEKYNLFE SSIRKTLPPA PEIKKRRPPP PIPPEENTE EIVVAMYDFQ ATEAHLRLLE RGQEYIILEK NDLHWWRRARD KYGSEGYIPS NYVTGKKSNN LDQYEWYCRN TNRSKAEQLL RTEDKEGGFM VRDSSQPGLY TVSLYTKFGG EGSSGFRHYH IKETATSPKK YYLAEKHAFG SIPEIIEYHK HNAAGLVTRL RYPVSTKGKN APTTAGFSYD KWEINPSELT FMRELGSGLF GVVRLGKWRA QYKVAIKAIR EGAMCEEDFI EEAKVMMKLT HPKLVQLYGV CTQQKPIYIV TEFMERGCLL NFLRQRQGHF SRDMLLSMCQ DVCEGMEYLE RNSFIHRDLA ARNCLVNEAG VVKVSDFGMA RYVLDDQYTS SSGAKFPVKW CPPEVFNYSR FSSKSDVWSF GVLMWEIFTE GRMPFEKNTN YEVVTMVTGR HRLHRPKLAT KYLYEVMRLC WQERPEGRPS LEDLLRTIDE LVECEETFGR</p>
Source	E.coli
Target Names	Tec
Protein Names	Recommended name: Tyrosine-protein kinase Tec EC= 2.7.10.2
Expression Region	1-630
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	<p>This protein belongs to the Tec family of non-receptor protein-tyrosine kinases containing a pleckstrin homology domain. Tec family kinases are involved in the intracellular signaling mechanisms of cytokine receptors, lymphocyte surface antigens, heterotrimeric G-protein coupled receptors, and integrin molecules. They are also key players in the regulation of the immune functions. Tec kinase is an integral component of T cell signaling and has a distinct role in T cell activation. This gene may be associated with myelodysplastic syndrome.</p>
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.