



Recombinant Arabidopsis thaliana F-box protein At1g47810 (At1g47810)

Product Code	CSB-BP680888DOA
Abbreviation	At1g47810
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4PT00
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MADTEKSLQS LDPIVDVLF EIFLNLPKAF LARFVCVSKL WAKIIRNQDF IRSF SFRSFR ENKQHRLLFA FKNQIKGYQE NWFYFSKSTH VSTPLYEPNS LVLDPQNIK FEALVGEEPP DFESSTVCHL KKMRYQNPSY VHGLISFLYG EEQIICNPSI GKSINLPTLG SSETIIGSFL GYDPIHAQYK VLCLNKSSLQ FCGHQVLTG AQNCSWRMIQ CTPHYPGTT SVCIDGVLYY SASRGFTMHE PLNLVRFDLR TESLEIASVF PEDFKSSVKP SLINYRGKVA VFSKYFCGDF GLWVLEDAKK QQWSKEQRIS IRRGVIGRLP QGLQILGTSD MGEVIFAPNY F
Source	Baculovirus
Target Names	At1g47810
Protein Names	Recommended name: F-box protein At1g47810
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.