



Recombinant Danio rerio Protein FAM110B (fam110b)

Product Code	CSB-EP680965DIL-B
Abbreviation	fam110b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4QRD7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MPTETLQADS KPVGQGGAFSA SAVPLRILNK GPDYFRRQTE PNPKRLSAVE RLEADKAKYV KSQEVINAKQ EPIKPQVLAK PPVCPVPLKR GGGVGSGGLKA SNNNAKSDTC ATTTKRENLN LEILKNILNS SSSSEGAAPKG TSHKVGARSW APHRSSSELS PGCRSFTDSL KVPPSTQGRH SPQGGNLNLS RRLLEERSCE VDGERSPLHA SHSSSDIRRM CNGKPLRAAR SSSSSAPPLP PKPSAKALAA CDNAPCSPER ETLTSTAAEL ELGSSVARRP SLHRKSDLS DRYARAGADV ERFFNYCGLD PEELESVGVE SFARANSDIV SLNFRSASMI SSDCEQSRHS NDDLTDDEEDD AGERVPYGIS AVERNARIK WLYSIKQARE SQKVSHV
Source	E.coli
Target Names	fam110b
Protein Names	Recommended name: Protein FAM110B
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.